



**ResMed**

**AirMini™**



Rewrite your  
bedtime stories



[ResMed.com/ExploreAirMini](https://www.ResMed.com/ExploreAirMini)



Actual



size

# What sets AirMini apart?

In addition to being the smallest CPAP on the market<sup>1</sup>, there are a number of reasons AirMini isn't just another travel CPAP.




## ActiveAir™

This crucial technology caps therapy flow while still increasing therapy pressure, which allows the system to be smaller, quieter and more efficient without compromising on therapy delivery or pressure swings.



## HumidX™

Our waterless humidification system allows exhaled air to pass through the HumidX before reaching the mask vent. The HumidX acts as a source of moisture throughout therapy, providing a comfortable portable therapy experience.<sup>2</sup>



## AirMini app

Patients can download the AirMini app on their smart device to control a range of comfort settings. Users receive their score each morning and can access the previous 30 days of sleep data, helping them stay engaged with therapy.

# The portable CPAP you

Finally, a portable CPAP you can trust! AirMini is the quality of therapy patients trust in our market. Leveraging the same proven algorithms and a rapid and effective therapy, AirMini sets a new bar for CPAP. 80% of users preferred AirMini over the two main



u've been waiting for

a complementary AutoSet™ device that delivers  
t-leading AirSense™ 10 bedside devices.

ange of new technologies that ensure discreet  
portable CPAP performance. In fact, more than  
rket-leading travel CPAPs in a recent study.<sup>3</sup>



# Optimized masks for AirMini

AirMini has been designed to work with our all-new AirFit™ 20 series and the AirFit P10 for AirMini masks. Patients can choose from full face, nasal and nasal pillows options – each designed to deliver comfortable and stable therapy.



Compatible masks include:



AirFit™ F20



AirFit™ N20



AirFit™ P10  
for AirMini

# A complete system

AirMini is more than just a portable CPAP, it's a turn-key solution for therapy on the go. The AirMini device and setup pack, optimized AirFit™ masks and AirMini app give patients what they need to bring high-quality therapy along no matter where life takes them. Optional accessories like the AirMini travel bag make it even easier for patients to pick up and go.



**1** As of March 20, 2017 **2** Currently not available when using full face masks **3** Results of independent, 3-night US study of 32 patients using ResMed AirMini, HDM Z1 Auto®, and Transcend Auto™ Mini CPAP Jan-Feb 2017.

Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA +1 858 836 5000 or 1 800 424 0737 (toll free). See ResMed.com for other ResMed locations worldwide. ActiveAir, AirFit, AirMini, AirSense, AutoSet and HumidX are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. © 2017 ResMed. 1019942/2 2017-04

ResMed.com