

**NE**  
gen. anytime. anywhere.®



The Inogen One® G3 System (IS-300):

- ✓ 4.8 pounds (includes single battery)
- ✓ Flow settings from 1-5
- ✓ One solution for oxygen at home, away, and for travel, 24/7

Call today! We'll help you determine which Inogen One System is right for you.

Join the tens of thousands of oxygen users who have reclaimed their freedom and independence. Once you try it, you'll wonder why you would use anything else!

Getting your INOGEN ONE oxygen is quick and easy. Call the number below to talk to one of our trusted HME Providers.



www.inogen.com

MKT-0389B

Reclaim Your Independence With The  
**INOGEN ONE**  
Family of Products



oxygen. anytime. anywhere.®

**INOGEN**  
oxy



The Inogen One® G4 System (IS-400):

- ✓ Only 2.8 pounds (includes single battery)
- ✓ Flow settings from 1-3
- ✓ Approximately half the size of the Inogen One G3
- ✓ One solution for oxygen at home, away, and for travel, 24/7



# INOGEN ONE<sup>G3</sup>



*“The day my Inogen Concentrator was delivered I got my life back. I can be out shopping for hours, travel, make long car trips. So many things are now possible.” - Carol Z.*

<b>Weight</b>	4.8 lbs. (includes single battery)
<b>Size</b>	Length: 8.75" Width: 3.0" Height: 8.25" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5
<b>Power Sensing</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 4.5 hours Recharge approx. 2-4 hours with AC or DC power Duration (double battery): Up to 9 hours Recharge approx. 4-8 hours with AC or DC power
<b>Noise Level</b>	39 dBA*
<b>Warranty</b>	3 year standard warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 use

\*At setting 2  
\*\* See Inogen One® Limited Warranty Statement for complete warranty information  
Additional accessories sold separately



## Freedom and Independence

Restoring your freedom and control of your day is a big part of what makes the Inogen One Systems so special. You can power an Inogen One System by plugging it into a wall outlet, DC outlet in your car, or for hours at a time with each rechargeable battery. So go ahead, jump in the car to run errands or take a weekend trip to see family. Come and go as you please without having to switch equipment. Or take that long trip you thought you'd never again have the chance to take. The Inogen One Systems are made to follow you anywhere, whether it's around town, on the road, or on a cruise. Where you go is now up to you with Oxygen. Anytime. Anywhere.®

### A single solution that can do it all!

The Inogen One System is one solution for stationary, portable, and travel use, day and night. Plus, it never has to be refilled.

# INOGEN ONE<sup>G4</sup>



<b>Weight</b>	2.8 lbs. (includes single battery)
<b>Size</b>	Length: 5.91" Width: 2.68" Height: 7.2"
<b>Oxygen Flow</b>	Pulse dose delivery system Three flow settings from 1-3
<b>Power Sensing</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power: for mobile use in car/airplane
<b>Battery</b>	Duration (single battery): Up to 2.7 hours Recharge approx. 3 hours with AC or DC power Duration (double battery): Up to 5 hours Recharge approx. 5 hours with AC or DC power
<b>Noise Level</b>	40 dBA*
<b>Warranty</b>	3 year standard warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 Use

\*At setting 2  
\*\* See Inogen One® Limited Warranty Statement for complete warranty information  
Additional accessories sold separately