

NE
gen. anytime. anywhere.®



The Inogen One® G5 System
(IS-500):

- ✓ 4.7 pounds (includes single battery)
- ✓ Flow settings from 1-6
- ✓ One solution for oxygen at home, away, and for travel, 24/7
- ✓ Inogen Connect Bluetooth App Included
- ✓ Meets FAA guidelines for use in flight

Call today! We'll help you determine which Inogen One System is right for you.

Join the tens of thousands of oxygen users who have reclaimed their freedom and independence. Once you try it, you'll wonder why you would use anything else!

Getting your INOGEN ONE is quick and easy.
Call the number below to talk to one of our trusted HME Providers.

Reclaim Your Independence With The
INOGEN ONE
Family of Products



oxygen. anytime. anywhere.®



inogen
www.inogen.com

MKT-0448A

INOGEN
oxy



The Inogen One® G4 System
(IS-400):

- ✓ Only 2.8 pounds (includes single battery)
- ✓ Flow settings from 1-3
- ✓ One solution for oxygen at home, away, and for travel, 24/7
- ✓ Inogen Connect Bluetooth App included
- ✓ Meets FAA guidelines for use in flight

INOGEN ONE^{G5}



"The day my Inogen Concentrator was delivered I got my life back. I can be out shopping for hours, travel, make long car trips. So many things are now possible." - Carol Z.

FREEDOM TO BREATHE



Freedom and Independence

Restoring your freedom and control of your day is a big part of what makes the Inogen One Systems so special. You can power an Inogen One System by plugging it into a wall outlet, DC outlet in your car, or for hours at a time with each rechargeable battery. So, go ahead, jump in the car to run errands or take a weekend trip to see family. Come and go as you please without having to switch equipment. Take that long trip you thought you'd never have the chance to take. The Inogen One Systems are made to follow you anywhere, whether it's around town, on the road, or on a cruise. Where you go is now up to you with Oxygen. Anytime. Anywhere.*

A single solution that can do it all!

The Inogen One System is one solution for stationary, portable, and travel use, day and night. Plus, it never has to be refilled.

Weight	4.7 lbs. (includes single battery)
Size	Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery)
Oxygen Flow	Pulse dose delivery system Six flow settings 1, 2, 3, 4, 5, 6
Power Sensing	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 6.5 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours Recharge up to 6 hours with AC or DC power
Noise Level	38 dBA*
Warranty	5 year standard warranty**
Operation	Simple control functions and easy-to-read LCD display Inogen Connect App
Use	Designed for 24/7 use

*At setting 2
Additional accessories sold separately

INOGEN ONE^{G4}



Weight	2.8 lbs. (includes single battery)
Size	Length: 5.91" Width: 2.68" Height: 7.2"
Oxygen Flow	Pulse dose delivery system Three flow settings from 1-3
Power Sensing	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power: for mobile use in car/airplane
Battery	Duration (single battery): Up to 2.7 hours Recharge approx. 3 hours with AC or DC power Duration (double battery): Up to 5 hours Recharge approx. 5 hours with AC or DC power
Noise Level	40 dBA*
Warranty	3 year standard warranty**
Operation	Simple control functions and easy-to-read LCD display Inogen Connect App
Use	Designed for 24/7 Use

*At setting 2
**See Inogen One® Limited Warranty Statement for complete warranty information