



# Golden Power Wheelchairs

HIGHLY DEPENDABLE WITH MAXIMUM MANEUVERABILITY,  
SUPERIOR STABILITY AND GREATEST FUNCTIONALITY

**WINTER 2018**



 GOLDEN



## THE MOST STABLE BASE IN ITS CLASS!

### YOU CAN'T STAND LIKE THIS ON ANY OTHER POWER WHEELCHAIR BASE

The Golden Compass™ line of power wheelchairs provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!

### THE NEW PORTABLE LITERIDER ENVY'S DESIGN MAKES OUR POWER WHEELCHAIRS THE EASIEST TO TRANSPORT!

**LiteRider**®  
GOLDEN  
ENVY  
2018

Available in 6 new vibrant  
colors: Cherry Red, Paradise Blue,  
Solar Flare Orange, Envy Green,  
Sunburst Yellow, Satin Silver



- Lightweight & Convenient
- Disassembles into 4 pieces
- Larger foot plate
- 22 Amp Hour battery for all-day range
- Off board battery charging for on the go



# Experience Constant Contact.

Inclines and thresholds are no match for the Golden Compass Sport's **true center wheel drive design with the patented springless articulating chassis!**

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving all Compass™ models unparalleled curb and threshold climbing ability. The Compass's suspension system pushes the center drive wheel down when the front and rear casters are raised and pushes the front and rear casters down when the drive wheels are raised. No other power base design can match the constant wheel to ground contact of any Golden Compass power base!

Scan QR code to view the Compass product video.



## LiteRider<sup>®</sup> GOLDEN ENVY 2018

**The Envy is available in 6 New Vibrant Colors!**

**Available Now:** Cherry Red & Paradise Blue

**Ships in 2018:** Solar Flare Orange, Envy Green, Sunburst Yellow & Satin Silver

**CHERRY RED**



LIMITED EDITION  
**SUNBURST YELLOW**



LIMITED EDITION  
**SOLAR FLARE ORANGE**



LIMITED EDITION  
**ENVY GREEN**



**PARADISE BLUE**



LIMITED EDITION  
**SATIN SILVER**



# Live the Golden Life

Test drive the industry's most comprehensive line of consumer power wheelchairs today!



Shown in Solar Flare Orange

**INTRODUCING**

## LiteRider<sup>GOLDEN</sup> ENVY<sub>2018</sub>

- Model GP162
- **Retail Cash Power Chair!**
- Larger Footplate
- New Shroud Design
- Wt. Cap.: 300 lbs.
- Turning Radius: 28.5"
- Batteries: 2-22AH
- Single Post Seat Mount
- **DISASSEMBLES FOR TRANSPORT**
- 6 NEW Vibrant Colors



## golden COMPASS<sup>TM</sup> SPORT

- Model GP605
- NEW aluminum wheels provide increased durability
- Medicare Codes: K0823 & K0822
- Wt. Cap.: 300 lbs.
- Turning Radius: 19.5"
- Batteries U1
- Single Post Seat Mount
- Dynamic Linx Controller Standard



<sup>GOLDEN</sup>  
*Alanté*  
**SPORT**

- Model GP205
- Medicare Codes: K0823 & K0822
- Wt. Cap.: 300 lbs.
- Turning Radius: 25"
- Batteries: 2-U1
- Single Post Seat Mount



*Alanté*  
BY GOLDEN

- Model GP207
- Medicare Codes: K0816 and K0815
- Wt. Cap.: 300 lbs.
- Turning Radius: 25"
- Batteries: 2 U-1
- Single Post Seat Mount



golden  
**COMPASS™**  
**HD**

- Model GP620
- Medicare Codes: K0825 & K0824
- Wt. Cap.: 450 lbs.
- Turning Radius: 20"
- Batteries NF22
- Four Post Seat Mount with Slide Seat Bracket



## POWER WHEELCHAIR PRODUCT SPECIFICATIONS

| Model Number                  | GP605<br>Compass Sport   | GP620<br>Compass HD               | GP162<br>LiteRider Envy                        | GP205<br>Alanté Sport       | GP207<br>Alanté             |
|-------------------------------|--------------------------|-----------------------------------|--|-----------------------------|-----------------------------|
| Medicare Code                 | K0823/K0822              | K0825/K0824                       | Not Coded                                      | K0823/K0822                 | K0816/K0815                 |
| Weight Capacity               | 300 lbs.                 | 450 lbs.                          | 300 lbs.                                       | 300 lbs.                    | 300 lbs.                    |
| Drive Wheels                  | Center                   | Center                            | Rear   | Front                       | Front                       |
| Maximum Speed                 | 4 mph                    | 4 mph                             | 3.5 mph  | 3.5 mph                     | 3.5 mph                     |
| Operating Range <sup>1</sup>  | 19.5 miles               | 12 Miles                          | 15.5 miles                                     | 16.5 miles                  | 16.5 miles                  |
| Ground Clearance <sup>3</sup> | 3.25"                    | 3.5"                              | 2.5"   | 1.75"                       | 1.75"                       |
| Turning Radius                | 19.5"                    | 20"                               | 28.5"  | 25"                         | 25"                         |
| Type Batteries                | 2-U1                     | 2-NF22                            | 2-22 AH  | 2-U1                        | 2-U1                        |
| Colors                        | Red, Blue                | Red, Blue                         | Red, Blue,<br>Orange, Yellow,<br>Green, Silver | Red                         | Red                         |
| Freewheel Mode                | Yes                      | Yes                               | Yes  | Yes                         | Yes                         |
| Electronic Speed Control      | Yes                      | Yes                               | Yes  | Yes                         | Yes                         |
| Electro-Mechanical Brakes     | Yes                      | Yes                               | Yes  | Yes                         | Yes                         |
| Charger                       | Off Board                | Off Board                         | Off Board                                      | Off Board                   | Off Board                   |
| Controller                    | Dynamic Linx<br>40 Amp   | Dynamic Shark<br>75 Amp<br>Remote | PG VR-2<br>50 Amp<br>Remote                    | PG VR-2<br>50 Amp<br>Remote | PG VR-2<br>50 Amp<br>Remote |
| Length                        | 39"                      | 40"                               | 36"  | 39.5"                       | 39.5"                       |
| Width                         | 24"                      | 25"                               | 22.5"  | 22.25"                      | 22.25"                      |
| WARRANTY                      | Frame                    | Lifetime                          | 5 Years  | 5 years                     | 5 years                     |
|                               | Drive Train              | 13 months                         | 13 months                                      | 13 months                   | 13 months                   |
|                               | Electronics <sup>2</sup> | 13 months                         | 13 months                                      | 13 months                   | 13 months                   |
|                               | Battery                  | 12 months                         | 12 months                                      | 12 months                   | 12 months                   |

1. Operating range will vary due to rider weight, drive surface and terrain.
2. Electronics warranty excludes batteries.
3. Measured at the lowest point mid-frame.

Note: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.