



Golden Scooters

SUPERIOR DESIGNS | FIRST CLASS LUXURY

WINTER 2018



Shroud Color Options Available



Silver



Pink



Yellow



BUZZY XLTM AROUND XLTM SERIES

COMPACT PORTABLE

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers more standard features to make your riding experience better than ever. Featuring a wrap around Delta Tiller, a new LED battery gauge, and a new ultra bright LED angle-adjustable headlight, the new Buzzaround XL is also longer to offer more foot room, with a great turning radius. The Buzzaround XL disassembles in seconds with just 4 easy steps making transporting it from one destination to another a breeze!



Optional color shrouds available.
*Upcharge applies.



GB107

BUZZY AROUND XLTM

BUZZAROUND XL

- 300 lb. weight capacity
- Increased foot room
- Color through chip resistant plastic (red only)
- Charger port in the tiller
- Quickly disassembles into five, lightweight pieces for easy transport
- Comes standard in matte red
- Pink, yellow and silver shrouds available
- Accepts new interchangeable color shrouds

**"Life is just better
on a Buzzaround!"**

(Buzz XL model shown below)





GB117D

**BUZZVI
AROUND XL™**

Optional color shrouds available.
*Upcharge applies.



Optional color shrouds available.
*Upcharge applies.



GB117H

**BUZZVI
AROUND XL
HD**

Optional color shrouds available.
*Upcharge applies.



**BUZZVI
AROUND XL™**



GB147D

Optional color shrouds available.
*Upcharge applies.



**BUZZVI
AROUND XL
HD**



GB147H

BUZZAROUND XL MODELS FEATURE:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport

WARRANTIES:

FRAME: Lifetime Warranty
DRIVE TRAIN: 2 Year Warranty
ELECTRONICS: 2 Year Warranty
BATTERIES: 1 Year Warranty



ALL BUZZAROUND XLHD MODELS FEATURE:

- Increased weight capacities up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red & Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22Amp batteries standard
- Operating range up to 14 miles per charge

WARRANTIES:

FRAME: Lifetime Warranty
DRIVE TRAIN: 2 Year Warranty
ELECTRONICS: 2 Year Warranty
BATTERIES: 1 Year Warranty



**BUZZYLS
AROUNDXLS™**
Comfort Spring Suspension



GB117S

**BUZZYLS
AROUNDXLS™
HD**



GB117Z

**BUZZYLS
AROUNDXLS™**
Comfort Spring Suspension



GB147S

**BUZZYLS
AROUNDXLS™
HD**



GB147Z

BUZZAROUND XLS MODELS FEATURE:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- NEW Comfort Spring Suspension for a smoother ride

WARRANTIES:

FRAME: Lifetime Warranty

DRIVE TRAIN: 2 Year Warranty

ELECTRONICS: 2 Year Warranty

BATTERIES: 1 Year Warranty



BUZZAROUND XLSHD MODELS FEATURE:

- Increased weight capacity up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22 Amp batteries standard
- Operating range up to 14 miles per charge
- NEW Comfort Spring Suspension for a smoother ride

WARRANTIES:

FRAME: Lifetime Warranty

DRIVE TRAIN: 2 Year Warranty

ELECTRONICS: 2 Year Warranty

BATTERIES: 1 Year Warranty



BUZZAROUND EXTREME

FULL-SIZE PORTABLE

The all new Golden Buzzaround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter. Both the three and four-wheel BuzzAround Extreme Scooters come standard with 2-U1 batteries – industry leading and the largest battery in this category – for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, and comfort suspension springs for the smoothest ride. Go the extra mile on a Buzzaround Extreme Scooter!

BUZZAROUND EXTREME FEATURES:

- 330 lb. weight capacity
- 4" ground clearance
- Full size & portable, yet easy to disassemble
- Larger 90 amp controller
- Larger 320 watt motor
- Larger 9" tires
- 2-U1 batteries for maximum range
- Ultrabright LED headlight and tail light
- Comfort suspension springs
- 18 mile range on one charge

**Largest batteries
and best range
in its class!*



GB118

WARRANTIES:
FRAME: Lifetime Warranty
DRIVE TRAIN: 2 Year Warranty
ELECTRONICS: 2 Year Warranty
BATTERIES: 1 Year Warranty



GB148



CHARGING DOCK:

The new Battery Docking Station for models GB118 and GB148 is portable and charges the interlocking U1 batteries. Perfect for keeping a set of charged batteries on hand, especially when visiting family or friends, or at a vacation destination.

BUMPER KIT

Can be ordered as a KIT or pre-installed on model GB118.



LiteRider^{GOLDEN}

LARGE PORTABLE

The unique tubular design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble, and maneuver smoothly indoors and out.



LiteRider^{GOLDEN}
LiteRider GL111D



LiteRider^{GOLDEN}
LiteRider GL141D

THE GL111D AND GL141D FEATURES:

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Adjustable tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat
- Adjustable armrests
- 2-22 Ah Batteries; increased operating range up to 12 miles
- Completely wireless disassembly
- 9" front and rear tires
- New Ultrabright LED headlight and rugged plastic basket standard!



LiteRider^{GOLDEN}

ENVY₂₀₁₈

LiteRider GP162 in Envy Green

Introducing the **NEW LiteRider Envy!** Available in 6 New Vibrant Colors with a Larger Foot Plate and New Shroud Design!

Golden's RETAIL CASH POWER WHEELCHAIR features:

- 300 lb. weight capacity
- 2-22 Ah Batteries; increased operating range up to 15.5 miles
- Turning radius: 28.5"
- Single post seat mount
- Disassembles easily for transport

CHERRY RED



PARADISE BLUE



LIMITED EDITION SOLAR FLARE ORANGE



LIMITED EDITION SUNBURST YELLOW



LIMITED EDITION ENVY GREEN



LIMITED EDITION SATIN SILVER



Red & Blue now available.
Other colors ship in 2018.

GOLDEN COMPANION

FULL-SIZE LUXURY

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:



Companion GC340
Shown with optional power lift seat.

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Backlit battery gauge
- Black, non-marking, low profile tires on sturdy steel rims
- Comfortable, fully adjustable padded armrests
- Under-carpet troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC240

Companion GC440

golden Avenger



Avenger GR541D

When you want a rugged, heavy duty scooters to tackle the challenging terrain of the great outdoors, look no further than the Golden Patriot™ and Avenger™. Both feature oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.

The Avenger model features:

- 500 lb. weight capacity
- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- 18-mile operating range

Golden Patriot™

The Patriot™ model features:

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge
- Black, pneumatic tires with brushed alloy wheels
- 13-Inch Pneumatic Tire and Rims



Patriot GR575D

AVAILABLE ACCESSORIES

*Only available on Companion models.



Oxygen Holder



Crutch Holder



Rear Basket



Safety Flag
(Included in Every Model!)



Power Elevating Seat*



Walker Holder



Cell Phone Holder




Cup Holder



Pack'n'Go Armbag



Pack'n'Go



Optional Buzzaround
Shroud Cover Colors.
Red & Blue come standard.



Quad Cane Holder



SPECIFICATIONS*

Item	Description	Length	Width	Front Tire	Rear Tire	Weight Capacity (lbs.)	Top Speed (mph)	Under Deck Clearance	Operating Range ¹ (miles)	Turning Radius	Batteries (2)	Weight w/ Batteries	Front Half Weight	Rear Half Weight	Seat Weight w/ Arms	Heaviest Piece	FreeWheel Mode	Electronic Speed Control	Charger
GB107D	Buzz LT (3-wheel)	40.5"	21.5"	8"	8"	300	4	3"	8'	37"	2-12Ah	98	33	25	21	33	Yes	Yes	Off-board
GB117D	Buzz XL (3-wheel)	42.5"	21.5"	8"	8"	300	4	3"	8'	37"	2-12Ah	101	35	25	20	35	Yes	Yes	Off-board
GB117S	Buzz XLS (3-wheel)	42.5"	21.5"	8"	8"	300	4	3"	8'	37"	2-12Ah	103	35	28	21	35	Yes	Yes	Off-board
GB117H	Buzz XLHD (3-wheel)	42.5"	21.5"	8"	8"	325	4	3"	14'	37"	2-22Ah	110	35	25	21	35	Yes	Yes	Off-board
GB117Z	Buzz XLSHD (3-wheel)	42.5"	21.5"	8"	8"	325	4	3"	14'	37"	2-22Ah	114	35	28	22	35	Yes	Yes	Off-board
GB147D	Buzz XL (4-wheel)	41.5"	21.5"	8"	8"	300	4	3"	8'	47"	2-12Ah	109	43	25	20	43	Yes	Yes	Off-board
GB147S	Buzz XLS (4-wheel)	41.5"	21.5"	8"	8"	300	4	3"	8'	47"	2-12Ah	111	43	28	21	43	Yes	Yes	Off-board
GB147H	Buzz XLHD (4-wheel)	41.5"	21.5"	8"	8"	325	4	3"	14'	47"	2-22Ah	118	43	25	21	43	Yes	Yes	Off-board
GB147Z	Buzz XLSHD (4-wheel)	41.5"	21.5"	8"	8"	325	4	3"	14'	47"	2-22Ah	122	43	28	22	43	Yes	Yes	Off-board
GB118	Buzz XLEX (3-wheel)	45"	22"	9"	9"	330	5	4"	18'	42"	2-U1	148	46	30	23	46	Yes	Yes	Off-board
GB148	Buzz XLEX (4-wheel)	44"	22"	9"	9"	330	5	4"	18'	56"	2-U1	155	53	30	23	53	Yes	Yes	Off-board
GL111D	LiteRider (3-wheel)	41.5"	22"	9"	9"	300	5	3.75"	12'	39"	2-22Ah	119	44	28	21	44	Yes	Yes	Off-board
GL141D	LiteRider (4-wheel)	43.25"	22"	9"	9"	300	5	3.75"	12.5'	48"	2-22Ah	127	52	28	21	52	Yes	Yes	Off-board
GP162	LiteRider Envy	36"	22.5"	N/A	9"	300	3.5	2.5"	15.5'	28.5"	2-22Ah	123	37	36	23	37	Yes	Yes	Off-board
GC240D	Companion (3-wheel)	42.5"	23.5"	9"	9"	350	4.5	3.75"	19.9'	39"	2-U1	166	46	48	28	48	Yes	Yes	Off-board
GC340D	Companion (3-wheel)	45.5"	23.5"	9"	10"	400	4.5	4.5"	20.9'	40.5"	2-U1	168	48	48	28	48	Yes	Yes	Off-board
GC440D	Companion (4-wheel)	47.5"	23.5"	10"	10"	400	4.5	4.5"	18.9'	53.5"	2-U1	183	63	48	28	63	Yes	Yes	Off-board
GA541D	Avenger	57.5"	24.5"	13"	13"	500	7	4"	18 ²	67"	2-NF22	265	68	68	55	68	Yes	Yes	Off-board
GR575D	Patriot	57"	28"	13"	13"	400	8	6.75"	24'	76.5"	2-75Ah	308	Base 152	N/A	50	160	Yes	Yes	Off-board

* ALL WEIGHT SHOWN IN LBS.

¹ Battery range tested at 200lbs., and will vary due to rider weight, drive surface and terrain.

² Battery range tested at maximum weight capacity, and will vary due to rider weight, drive surface and terrain.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems. **Rear wheel drive:** All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

Choose from a Complete Line of Lifestyle Solution Products by Golden



POWER WHEEL CHAIRS

- Compass Sport
- Compass HD
- LiteRider Envy
- Alanté
- Alanté Sport



POWERED LIFT AND RECLINE CHAIRS

- Value Series
- Signature Series
- Comforter Series
- Comforter Wide Series
- MaxiComforter Series



Visit Your Local Authorized Golden Dealer.



www.goldentech.com