


INOGEN ONE^{G5™} System



Inogen One
AC Power Supply*
Item #BA-501



Inogen One G5
DC Power Cable*
Item #BA-306



Inogen One G5
Lithium Ion Battery*
Up to 6.5 hours run time
Item #BA-500



Inogen One G5
Carry Bag*
Item #CA-500



Inogen One G5
Backpack**
Item #CA-550



Inogen One[®]G5
Extended Life Lithium Ion Battery**
Up to 13 hours run time
Item #BA-516

The Inogen One G5 System (IS-500):

- Compact and Lightweight
- Flow Settings 1-6
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use
- Easy-to-read LCD display
- Inogen Connect App

*Included with system

**Accessory sold separately

INOGEN ONE^{G5}

Oxygen. Anytime. Anywhere[®]

Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G5 offers the most oxygen per pound for a portable oxygen concentrator on the market today. With flow settings from 1-6, the Inogen One G5 is designed to dramatically increase independence for most supplemental oxygen users 24/7. Advanced diagnostics take the guesswork out of patient support, as the Inogen One G5 is bluetooth enabled and compatible with Inogen Connect, the free mobile app. It can be charged at home or on the go - giving patients the freedom of Oxygen. Anytime. Anywhere.

Weight	4.7 lbs. (includes single battery)
Size	Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery)
Oxygen Flow	Pulse dose delivery system Six flow settings 1, 2, 3, 4, 5, 6
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 6.5 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours Recharge up to 6 hours with AC or DC power
Noise Level	38 dBA*
Warranty	5 year**
Operation	Simple control functions and easy-to-read LCD display Inogen Connect App
Use	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

*At flow setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive
Goleta, California 93117
855-631-2438

www.inogen.com/hme-providers



MKT-0447A