

INDEPENDENCE™

Make the Best of Yours



PVI OUTSIDE VEHICLE LIFTS






Optional Battery Kit
Available for Every Lift

**Verify Vehicle
Compatibility at
www.pviramps.com**



Adjustable
Cradle.



INDE1

INDE2

Mobility Devices

Most 3 and 4 wheel scooters with a wheelbase up to 42"

Most Power Wheelchairs

Weight Capacity

350 lbs.*

350 lbs.*

Platform Size

27.5" x 39" w/ adjustable rear wheel cradle it is extended to 46"

27.5" x 39"

Installed Weight

95 lbs.

98 lbs.

Securement Method

Automatic hold down foot. As the user raises the platform up, the hold down foot will automatically come down and secure the scooter.

4 Q'Straint tie downs. Most dependable securement option, however, not ideal for those with limited dexterity.

Hitch Class

Class II or III

Class II or III

Benefits

Each Lift Features: - One Switch Operation - No Chair or Scooter

*Optional weight capacity upgrade to 400 lbs. available (Model # AL902).

Our Lifts Are Proudly Made In The USA.



Securement Methods



Automatic Hold-down Foot



Automatic Hold-down Arm



Q'Straint Retractable Hooks



INDE3

Mid-Wheel drive power wheelchairs with 10" drive wheels


350 lbs.*

16.5" x 27"

63 lbs.

Padded hold down arm. Provides quick and easy securement. Great for those with limited dexterity or other issues.

Class II or III



INDE4

Larger Power Wheelchairs

350 lbs.*

32.5" x 38.5"

115 lbs.

Padded hold down arm. Provides quick and easy securement. Great for those with limited dexterity or other issues.

Class III or greater

Modifications - Included License Plate Mount - 3-Year Warranty

INDEPENDENCE™ Accessories



Swing-Away - AL105



Battery Kit - ALBAT



Wheelchair and Scooter Covers

Also Available
from **PVI**



Portable Ramps



Threshold Ramps



Modular Ramps

	Model #	Description	Weight	MSRP		
Lifts	INDE1	Universal Scooter Lift	95 lbs.	\$1,999.00		
	INDE2	Universal Power Chair Lift	98 lbs.	\$2,199.00		
	INDE3	Compact Auto Power Wheelchair Lift	83 lbs.	\$2,299.00		
	INDE4	XL Auto Universal Power Chair Lift	115 lbs.	\$2,549.00		
Accessories	AL105	Swing-Away	43 lbs.	\$449.00		
	AL902	Weight Capacity Upgrade to 400 lbs.	N/A	\$199.00		
	ALHA	Hitch Adapter Universal	7 lbs.	\$179.00		
	AL4Q	4-Q'Straints	N/A	\$399.00		
	Covers	AL010C	Micro Scooter Cover	3 lbs.	\$149.00	
		AL030C	Universal Manual Wheelchair Cover	3 lbs.	\$149.00	
		AL050C	Micro Power Chair Cover	3 lbs.	\$149.00	
		AL100C	Universal Scooter Cover	4 lbs.	\$149.00	
		AL100CE	Value Line Scooter Cover	4 lbs.	\$99.00	
		AL500C	Universal Power Chair Cover	4 lbs.	\$149.00	
		AL500CE	Value Line Power Chair Cover	4 lbs.	\$99.00	
		AL500CXL	XL Power Chair Cover	6 lbs.	\$179.00	
	AL100CXL	XL Scooter Cover	6 lbs.	\$179.00		
Batteries	ALBAT	Battery Kit w/ A/C Adapter	2 lbs.	\$268.50		
	18VB	18V Lithium Battery (Replacement/Spare)	1.8 lbs.	\$328.90		
	110VC	110V Battery Charger (Replacement/Spare)	3.3 lbs.	\$164.90		
	110/12VDC	Dual Power Home & Vehicle Battery Charger	1.6 lbs.	\$286.90		

Making Access a Way of Life

2620 Industrial Drive, Fairbury, NE 68352 | 800-554-7267 | www.pviramps.com

