



PHILIPS

RESPIRONICS

SimplyGo

Life made
simple

Philips Respironics SimplyGo portable oxygen concentrator –
simply changes your patient's life

SimplyGo




The freedom your patients deserve.

A short shopping spree, a museum visit, a pleasant weekend, or a well deserved holiday— with SimplyGo, your oxygen-therapy patients may enjoy it all. Lightweight and durable, delivered with a handy carrying case and a cart with easy-to-go wheels, SimplyGo is always ready for a stroll. Its accessory case keeps everything on hand, for all occasions.



SimplyGo is Philips' first oxygen concentrator that offers the benefits of a continuous flow mode in a lightweight, easy-to-carry design of 4.5 kg. No need to compromise between an effective oxygen therapy output for your patients and the usual weight of the device. SimplyGo offers a solution that fulfills the needs of most of your oxygen-therapy patients. With SimplyGo, patients have the flexibility to be away for the day, for work, for a weekend or a long holiday.

Three oxygen-delivery modes are available to fit most patients' needs:

-  **Pulse mode**
SimplyGo detects breath intake and delivers the necessary volume of oxygen.
-  **Continuous flow**
Continuous oxygen delivery from 0.5 to 2 LPM.
-  **Sleep mode**
Pulse mode has an increased trigger sensitivity for use during sleep. SimplyGo switches automatically to 2 LPM *Continuous flow* if no breath intake is detected.

With oxygen always on hand, your patients have the freedom to go anywhere, anytime.

At home, the **EverFlo** stationary concentrator is the ideal solution to complement the **SimplyGo** ambulatory device for patients requiring home oxygen therapy. Since **SimplyGo** and **EverFlo** use **non-delivery technology**, they eliminate the inconvenience of recurrent oxygen deliveries and can help improve patients' independence and quality of life.*



*Portable oxygen increases the daily duration of oxygen therapy in patients with COPD and improves the quality of life of the most compliant patients. Vergeret et al. Eur Respir J 1989;2: 20-25

SimplyGo

Reliability and maintenance made simple!

COPD (chronic obstructive pulmonary disease) patients are dependent on oxygen therapy equipment that they need for many hours* a day. **SimplyGo** offers the independence and freedom patients want in a durable design that can resist the bumps and bruises of daily activities while providing an effective therapy.



Built-in reliability

SimplyGo strives to set a new standard in portable concentrators' reliability. **SimplyGo** has been designed, built, and rigorously tested to withstand the stress of everyday life when following patients anytime, anywhere.

During development, SimplyGo was subjected to extreme conditions including impact, vibrations, and temperature extremes to test reliable performance day-in and day-out. With 15000 hours of compressor lifetime, its performance is trustworthy. This may mean fewer callbacks for you and fewer frustrations for your patients due to accidental damage.

You can count on SimplyGo.

Maintenance made simple!

SimplyGo makes it easy for providers, as it requires virtually no preventive or patient maintenance.

SimplyGo is flexible:

- Make a choice between two-, three-, or five-year warranty. Pick the plan that's best for your service.
- SimplyGo's easily accessible parts and comprehensive service tools allow your technicians to service the device in-house instead of returning it to the manufacturer.

Notes

*Clinical guidelines recommend that COPD patients should use oxygen therapy for at least 15 hours per day to have positive effects on their long-term survival and their quality of life. Therefore, these patients require oxygen at home (stationary oxygen therapy) and outside (ambulatory oxygen therapy).



SimplyGo

At-a-glance



Simple-to-select oxygen delivery mode

Simple-to-use touch-pad interface

Lightweight - 4.5 kg
The only lightweight device able to deliver both, pulse and continuous oxygen flow.

Simple-to-carry anywhere, anytime!

Simple-to-replace/remove rechargeable battery.
See product specifications table for more information about battery life.

The United States Federal Aviation Administration (FAA) has approved the SimplyGo for in-flight use on commercial airliners.¹



SimplyGo

Ready for a stroll



SimplyGo sturdy cart with large wheels for easy stroll

Humidifier pouch

DC power supply

AC power supply

External battery module

Airline power supply

External battery charger

All your patient needs in one SimplyGo accessory bag

¹Although the FAA and other international aviation regulatory bodies allow the use of approved portable oxygen concentrators on commercial airlines, many domestic and international carriers have additional requirements, including advance notification, and possession of a physician's letter and oxygen prescription.

SimplyGo specifications

Product specifications

Oxygen concentration ²	87 to 96% at all settings	
Flow settings		
Pulse-mode bolus size	1 = 12 ml 2 = 24 ml 3 = 36 ml	4 = 48 ml 5 = 60 ml 6 = 72 ml
Continuous-flow mode	0.5 LPM - 1 LPM - 1.5 LPM - 2 LPM	
Night mode	Variable pulse volumes based on breath rate, to maintain a constant minute volume per setting: 1 = 250 ml/min 2 = 500 ml/min 3 = 750 ml/min 4 = 1000 ml/min 5 = 1250 ml/min 6 = 1500 ml/min	
Sound level	43dB (nominal) at pulse setting 2	
Alarms	Low oxygen purity (<83%) No breath High breath rate No flow Low battery Depleted battery System technical fault	
Dimensions	29.2 cm x 25.4 cm x 15.2 cm 11.5 in x 10 in x 6 in	
Weight	4.5 kg 10 lbs with battery installed	
Power consumption	150W while charging	120W while not charging

Operating conditions

Operational temperature	5° to 40°C 41° to 104°F
Relative humidity	15 to 95%
Altitude	Up to 3048 m 10000 ft

Battery autonomy

Pulse mode setting 2*	3 hr min
Continuous flow 2 LPM	54 min

*Based on 20 breaths per minute

²Based on an atmospheric pressure of 14.7 PSIA (101 kPa) at 21°C | 70°F

Ordering information

SimplyGo international system

	Part number
Includes: SimplyGo portable oxygen concentrator, carrying case, mobile cart, one rechargeable Lithium ion battery, DC and AC power cords, accessory case and user manual	1069058

Accessories	Part number
Carry cart	1074885
SimplyGo carrying case	1082663
Lithium ion battery	1082662
Power cord (Standard Europe)	1082665
External battery charger (European plug)	900-118
External battery charger (UK/Ireland plug)	900-119
Accessory case	1083696
AC supply	1082661
DC car charger	1083692
DC airline cord	1083693
Humidifier pouch	1083699
Humidifier connecting tube	1039642
External battery module	1109643
Liberty kit (external battery module and SimplyGo battery)	1116269

Warranty options	Part number
SimplyGo two-year total warranty	1097623
SimplyGo three-year total warranty	1097577
SimplyGo five-year total warranty	1097778



How to reach us:
www.philips.com/healthcare
healthcare@philips.com

Product information:
www.philips.com/simplygo

Philips Respironics, SimplyGo, and EverFlo are trademarks of Koninklijke Philips N.V. and its affiliates.
All rights reserved.

©2014 Koninklijke Philips N.V. All rights reserved.
Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.